

Introducing X-halo

X-halo is the first clinically proven device which allows asthmatics to track airway inflammation at home.

X-halo is quick and easy. It takes less than 2 minutes a day to use and is suitable for everyone, children and seniors alike.

Airway Inflammation & Asthma

Airway inflammation occurs in all asthma attacks but did you know that in 3 out of 4 cases*, the inflammation can begin to develop unnoticed up to 48 hours before the actual attack?


The revolutionary new X-halo allows asthmatics to track their airway inflammation levels every day, at home, in just minutes.



**Caused by bacterial or viral infections*

Airway Inflammation & X-halo

In its early stages, the inflammation is undetectable to asthmatics, leaving them unaware of the growing threat.



As a result, they are not given the opportunity to take preventive action.

Until now.

X-halo allows asthmatics to track their airway inflammation levels easily and quickly at home. Giving them the opportunity to take more control over their asthma management.



How Does X-halo Work?

X-halo precisely measures the temperature of the air in your lungs - called Exhaled Breath Temperature (EBT) - to within a tenth of a degree. When your lungs begin to get inflamed, your EBT rises and X-halo lets you know.

X-halo[™]
Breathe easier



DMi
DELMEDICA

Dedicated to improving lives of people
with meaningful innovations in respiratory care

www.x-halo.com

02022015

The X-halo App: Complete Asthma Management

X-halo's custom-designed application acts as the brain of X-halo and as it gets to know the user over time, it can be used to identify trends and patterns, empowering asthmatics to make more informed and proactive asthma management choices.

It even goes beyond airway inflammation tracking to integrate all aspects of asthma management into the user's daily routine:

- Asthma Action Plan upload
- Monthly Asthma Questionnaire
- Medication tracking



X-halo Improves The Lives of Asthmatics

With X-halo, asthmatics finally have a way to track the state of their lungs on a daily basis, at home.

When used regularly X-halo can improve the lives of asthmatics in many ways:

- More informed lifestyle choices
- Peace of mind
- Increased asthma control
- Better asthma management



Power supply	2 AAA batteries
Battery life	120 measurements or 2 months under normal use
Connections	Connects to smart phone, tablet, PC or Mac through an audio jack
App	Available free from the App store for iOS devices or Google play for Android devices and downloadable from www.x-halo.com for your PC or Mac
Measurement	Range: 25-42°C
Accuracy	Sensor accurate to 0.1°C in the temperature range specified
Storage	Temperature range: -20-60°C; Humidity range: 10-90%RH
Usage	Temperature range: 10-35°C; Humidity: 30-80% RH
Characteristics	Non-predictive, no measurement offset and no-calibration necessary

X-halo™



A revolution in home asthma management